



Vis Pro

Ophthalmic Surgical Operating Microscope

Sanma Vis Pro — Designed for the Depths of
Retinal Precision



Deep Focus,
Wide **Vision**



Ultra-Bright Posterior Illumination

Uniform illumination engineered to maintain clarity during vitrectomy, endo-laser, membrane peeling, and fluid-air exchange. Adjustable LED illumination minimises glare and retinal phototoxicity.

Widefield Fundus Optics

High-NA objective designed for full posterior visualization, perfect with BIOM-type add-on systems.

Deep Depth of Field

Stable focus across shifting retinal planes, reducing manual focusing effort during dynamic surgical steps.

Motorised XY control and focus enable precise micro-movements without repositioning

Sanma Vis Pro — Designed for the Depths of Retinal Precision

Apochromatic Clarity in Every Layer

- True-color retinal imaging
- High contrast for ILM/ERM differentiation
- Anti-reflective coatings for glare-free views

Precision at Your Foot Pedal

- Smooth motorized zoom & micro-focus
- Seamless switching between widefield and high magnification
- Intuitive tactile response designed for retinal workflow

Documentation & Teaching Excellence

- Integrated HD camera
- Real-time recording & assistant view
- Ideal for teaching hospitals & complex surgical demonstrations

Ergonomics for Long Retinal Cases


- Balanced arm system
- Long-range articulation
- Reduces surgeon fatigue over extended procedures

Designed for the Depths of Retinal Precision

Engineered for Depth, Detail & Delicate Retinal Work

Paragraph in Graphite Grey

Sanma VisPro delivers the high-contrast, wide-field, and deep-focus visualisation required for vitreoretinal surgery. Designed for prolonged procedures and fine-scale retinal manipulation, the system provides unwavering clarity, surgeon comfort, and precise illumination within the posterior segment.



VisPro

Features	Description	Configurable
Optics	Fully apochromatic Galilean type, broad bandwidth anti-reflective coated	
Binocular tube	Binocular tube: 45° inclined binocular tube, focal length f=170mm, IPD adjustment of 52-76mm.	Binocular tube: 0° - 210° tiltable binocular tube, focal length f=170mm, IPD adjustment of 52-76mm.
Stereo base	22mm	
Eye piece	10x pair of wide field push sleeves & high eye point	
IPD	Inter pupillary distance: 50mm to 75mm	
Diopter adjustment	-6D to +6D	-5D to +8D
Magnification	Motorised zoom 1:6 ratio	
Objective Working Distance	Apochromatic Optics 200mm.	
Motorised focusing	Min: 30mm / Max: 50mm	
Illumination	Illumination: Coaxial LED.	
Field of Illumination	Ø70.5mm	
Intensity control	Continuously adjustable lamp intensity control in handgrip buttons	
XY	X-Y coupling range minimum 40x40mm with XY reset button.	
Foot Pedal	Motorised foot switch control for XY, Zoom, focus, illumination and red-iglow on/off, recorder on/off. (water proof).	
Arm Movement	Articulated spring balance, well counter balance for smooth positioning	
Floor Stand	Stable and sturdy rust free floor stand on 4 smooth castor wheels with 2 lockable	
Imaging System	<ul style="list-style-type: none"> Integrated full HD Camera with Video recording & still image storage 7" TFT Monitor integrated with carrier arm 	
Power	AC 100 – 120V/220-240V, 50-60Hz 280V	
Asepsis	Sterilisable autoclave able silicon caps for all control knob - 2 Sets.	

Configurable Imaging System



FULL HD
1928x1080

4K
3840 x 2160



Live Video &
Still image



Monitor
27" / 32"



Sanma
VisPro

Compatibility with Accessories

- Works seamlessly with BIOM-type wide-angle viewing systems
- Compatible with endoilluminator and 3D visualisation add-ons
- 45° or 0° - 210° tilt able binocular tube

Sanma's 3D Camera System in Cataract & Retinal Surgical Procedures

Enhanced Depth Perception

- Provides optimal natural stereoscopic depth.
- Surgeons can assess spatial relationships of ocular structures more accurately.
- Essential for precision in capsulorhexis, nucleus manipulation, membrane peeling, and vitreous base shaving.

Superior Teaching & Training Experience

- Entire OR team sees the same 3D view as the surgeon.
- Trainees understand instrument handling, tissue planes, and surgical flow more intuitively.
- Allows teaching of complex steps in retinal surgery (e.g., ILM peeling) with unmatched clarity.

Ergonomically Comfortable Surgery

- Surgeons operate looking straight at a 3D monitor, reducing neck and back strain.
- Longer retinal procedures become more comfortable and stable.
- Minimizes fatigue, supporting higher surgical precision.

Digital Image Enhancement in Real Time

- Adjustable contrast, brightness, color saturation, and digital zoom.
- Enhances visibility in low-light or media-opacified situations.
- Useful in retinal surgeries requiring delicate visualization—especially macular work and membrane differentiation.

Wider Field of View for the Surgical Team

- Everyone—assistant, scrub nurse, anesthetist—shares the same high-definition 3D picture.
- Improved coordination during critical steps such as IOL insertion, endolaser, or fluid-air exchange.

Improves Documentation & Case Review

- High-resolution 3D video recording.
- Enables post-operative review, academic presentations, and digital archiving.
- Ideal for teaching hospitals and surgical training programs.



Sanma 3D Visualisation Elevates Every Eye Surgi Procedure



Specifications

- Sensor: 2X 1/1.3" CMOS sensor
- Resolution: 4K UHD - 3840x2160
- Video output: 2X 3G-SDI
- 3D output: side by side
- Recording format: 3D & 2D, MP4 / MOV
- Still image: 2D / 3D .jpeg, PNG
- HDD: 4TB
- Motorised adjustable monitor stand

■ Monitor* : 32" / 55"



Sanma VisPro

Ophthalmic Surgical Operating Microscope

EMC 2014/108/EC
EAR Regulation (EU) 2017/745 on Medical Devices
Clinical Evaluation report as per EU MDR 2017/745
Technical documentation as per Annex II of EU regulation 2017/745
ISO 13485:2016
CDSCO: TN/M/MD/OO4755


Specifications & Design are subject to change without prior notice as a result of further product development
©Sanma Medineers Vision Private Limited
All Rights Reserved
2025



www.sanmas.com

Sanma Vis Pro – Retina Surgery Excellence



 Sanma Medineers Vision PVT LTD
68-A, Aavin Dairy Road, North Phase,
SIDCO Industrial Estate, Ambattur,
Chennai, TN - 600098, India

 Email: sales@sanmas.com